

CeloNova introduces

# Embozene TANDEM™

(Drug-Elutable Microspheres)

Targeted Embolization Made EZ

- TANDEM™ combines the precision and ease of use you have always appreciated in Embozene® Microspheres **with full integration of drugs.**\*
- Capable of loading doxorubicin-HCl and irinotecan-HCl up to 50 mg/ml Microspheres.

Efficiency | Precision | Flexibility

# TANDEM™





**For health care providers, some benefits of choosing TANDEM™ include**

**Efficiency**

- Load up to 100 mg of drugs in a 2 ml syringe
- Load up to 150 mg of drugs in a 3 ml syringe
- Fast drug loading time:
  - Doxorubicin: 60 minutes
  - Irinotecan: 30 minutes
- Controlled drug release rates for optimal localized drug delivery

**Precision**

- Targeted drug delivery with tightly calibrated Microspheres and superselective TACE
- Controlled and precise drug delivery and release
- Microsphere sizes remain stable during drug loading and storage. Typical size change is less than 5%

**Flexibility**

- Ideal for passage through microcatheters

**Due to the design characteristics Embozene TANDEM™ Microspheres may be loaded with\***

- Doxorubicin-HCl for local, controlled, sustained dose elution to targeted tumor sites such as hepatocellular carcinoma (HCC) after embolization
- Irinotecan-HCl for local, controlled, sustained dose elution to targeted tumor sites such as metastatic colorectal cancer (mCRC) after embolization

**TANDEM™ Microspheres are available in the following sizes**

Size	2 ml Syringe	3 ml Syringe
40 ± 10 µm	10420-TS0	10430-TS0
75 ± 15 µm	10720-TS0	10730-TS0
100 ± 25 µm	11020-TS0	11030-TS0

To order Embozene TANDEM™ please contact your local CeloNova representative.

[www.celonova.com](http://www.celonova.com)

\*Embozene TANDEM Microspheres are indicated for embolization according to the IFU. They may be loaded with a drug to elute a local, controlled, sustained dose of said drug to targeted tumor sites after embolization. Loading of drug should be under a physician's direction, choice and responsibility, based on type and dose of drug most beneficial to the patient.

\*This product is currently available outside the U.S. only.

The CeloNova BioSciences, and Embozene TANDEM™ names and logos are registered to CeloNova BioSciences. All product names in any character format and/or with the trademark symbol are trademarks of CeloNova BioSciences, Inc. Manufactured at CeloNova BioSciences Germany GmbH, Ulm, Germany (CE0297). Distributed by CeloNova BioSciences, Inc. San Antonio, USA (CE0086). ©2012 CeloNova BioSciences. All Rights Reserved.



CE0086 BIOSCIENCES

APM0016-06 130115 - Rev A